



What's new in ChairsideCAD 3.0 Galway SR1

Release presentation

12/2021

exocad



Highlights for ChairsideCAD SR1

Mock-up tooth setups for improved patient communication

Intro


Design a mock-up tooth setup based on the pre-op model and print the mock-up model, or the mock-up itself, to create a clip-on custom try-in. This creates new possibilities for patient consultation with *Smile Creator*. Patients can now physically try out new smiles with their try-ins.

Technical description

- The additional function is available as part of the *ChairsideCAD 3.0 Galway Core* version with *Smile Creator*
- Users with ChairsideCAD versions activated before January 1, 2022 need to activate the free-of-charge *MockupTooth* module to create mock-ups

Benefits

- Create a clip-on mock-up of the planned anatomic situation in the workflow for a printable try-in that wraps around the patient's teeth
- Place finished product on patient's teeth for an initial fitting of the tooth shapes
- Possibility to print mock-up models / digital waxup models



Save insertion direction to scanning profile

Intro


Save the insertion direction of the restoration during the scan process and import it directly into the *ChairsideCAD* software.

Technical description

- Insertion direction is automatically applied in the design
- Additional information is saved in a scaninfo file
- Design step is no longer shown in the wizard


Benefit

- Dentists can fix insertion direction and margin line immediately after scanning and before transferring the file to the CAD software
- Fewer steps in the design software required



Additional languages

- Portuguese
- Spanish
- French





exocad

Thank You

exocad GmbH

[exocad.com](https://www.exocad.com)

 facebook.com/exocad

 facebook.com/groups/exocadexperts

 linkedin.com/company/exocad-gmbh

 videos.exocad.com